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### IgG Food Intolerance Test

Many people exhibit chronic food sensitivity. Unlike the immediate effects of IgE-mediated food allergy, IgG-mediated food intolerance reactions may take several days to appear. Symptoms include lethargy, weight gain, dermatitis, arthritis, migraines and irritable bowel syndrome (IBS). Controlled removal of the suspect foods from the patient's diet can, in many cases, rapidly improve their condition.

### FoodPrint® Technology

The FoodPrint® microarray is produced by 'printing' over 220 food extracts, in duplicate, onto nitrocellulose pads attached to a glass microscope slide, together with calibration standards and controls. Each protein spot measures only 130µm in diameter. Blood samples are diluted and dispensed onto the printed microarray and any food IgG antibodies present in the sample will bind to the food extracts. Bound food IgG antibodies are detected through the use of immunoassay reagents which produce a blue reaction in the presence of an antibody. The intensity of the reaction is measured using a high-resolution optical scanner and calibrated against known standards, using proprietary software to produce quantitative results. Concise, colour-coded personalised food reports are generated, allowing the patient or health professional to determine the most reactive foods with ease.

The FoodPrint® system has been installed in almost one hundred laboratories globally to-date. Hundreds of thousands of tests have been performed since its launch in 2007, providing patients with valuable dietary information that has enabled them to improve their health.

### Features & Benefits

- Works with whole blood, plasma or serum
- 5µl sample volume - finger-prick sample is less invasive for the patient
- High throughput multiplex assay - 64 samples assayed in 2 hours
- Highly miniaturised – 1 slide instead of 80 microtitre plates
- Cost savings – no need for microtitre plates and associated equipment
- Fully quantitative results – reactive foods can be easily identified
- Flexible array design allowing for dietary variation
- Detects over 220 food IgG antibodies with no pooled/mixed testing
- Wide range of food panels available – accommodates any budget or dietary need


### Assay Methodology

The FoodPrint® assay is simple to perform and all incubations are conducted at room temperature. It takes less than 2 hours to complete 1-4 slides containing 16-64 patient samples.

Cambridge Nutritional Sciences Ltd (CNS) has developed FoodPrint®, a laboratory-based system utilising state-of-the-art microarray technology. Originally invented for studying DNA and gene expression, CNS has extended the microarray technology to Food Intolerance testing. FoodPrint® is a sensitive and accurate technique that can detect the presence of IgG food-specific antibodies to over 220 commonly eaten foods.

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## Results

Fully quantitative results are generated for every patient tested. Results are presented in personalised food reports which identify the precise IgG antibody level for every food item selected. The following reports are also available from the FoodPrint® 200+ test: FoodPrint 40+, FoodPrint 60+, FoodPrint 120+, FoodPrint Vegetarian, FoodPrint Vegan and FoodPrint Herbs/Spices. Reports are concise, colour-coded and grouped according to food reactivity or food groups, allowing health professionals to easily identify 'problem' foods:

### Report 1

Foods listed alphabetically in food groups

### Report 2

Foods listed in order of reactivity

### Ordering Information

PRODUCT CODE	DESCRIPTION	PACK SIZE
CNS350	FoodPrint® Microarray System Package	1 system
CNS301-1/4	FoodPrint® Microarray 200+ Food IgG (4 patient)	4 tests
CNS301-4/4	FoodPrint® Microarray 200+ Food IgG (16 patient)	16 tests
CNS301-16/4	FoodPrint® Microarray 200+ Food IgG (64 patient)	64 tests

Kits to be stored at 2°C to 8°C. Shelf life - 18 months

All products are CE-marked against the relevant directive where required

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**FOODPRINT**

<b>Dairy/Egg</b> 	Alpha-Lactalbumin Beta-Lactoglobulin	Buffalo Milk Casein	Cow's Milk Egg White	Egg Yolk Goat's Milk	Sheep's Milk
<b>Fish &amp; Seafood</b> 	Anchovy Barnacle Bass Carp Caviar Clam Cockle Cod	Crab Cuttlefish Eel Haddock Hake Herring Lobster Mackerel	Monkfish Mussel Octopus Oyster Perch Pike Plaice Razor Clam	Salmon Sardine Scallop Sea Bream Shrimp/Prawn Sole Squid	Swordfish Trout Tuna Turbot Winkle
<b>Fruit</b> 	Apple Apricot Avocado Banana Blackberry Blackcurrant Blueberry Cherry	Cranberry Date Fig Grape Grapefruit Guava Kiwi Lemon	Lime Lychee Mango Melon Mulberry Nectarine Olive Orange	Papaya Peach Pear Pineapple Plum Pomegranate Raisin Raspberry	Redcurrant Rhubarb Strawberry Tangerine Watermelon
<b>Grains</b> 	Amaranth Barley Buckwheat Corn (Maize)	Couscous Durum Wheat Gliadin Malt	Millet Oat Polenta Rice	Rye Flour Spelt Tapioca	Wheat Wheat Bran
<b>Herbs &amp; Spices</b> 	Aniseed Basil Bayleaf Caromile Cayenne Cinnamon Clove	Coriander Cumin Curry Dill Garlic Ginko Ginger	Ginseng Hops Licuorice Marjoram Mint Mustard Seed Nettle	Nutmeg Parsley Peppercorns (B/W) Peppermint Red Chili Rosemary Saffron	Sage Tarragon Thyme Vanilla
<b>Meat</b> 	Beef Chicken Duck	Horse Lamb Ostrich	Ox Partridge Pork	Quail Rabbit Turkey	Veal Venison Wild Boar
<b>Nuts/Seeds</b> 	Almond Brazil Nut Cashew Nut	Coconut Flax Seed Hazelnut	Macadamia Nut Peanut Pine Nut	Pistachio Rapeseed Sesame Seed	Sunflower Seed Tiger Nut Walnut
<b>Vegetables</b> 	Artichoke Asparagus Aubergine Bean (Broad) Bean (Green) Bean (Red Kidney) Bean (White Haricot) Beetroot Broccoli	Brussel Sprout Cabbage Cabbage (Red) Caper Carrot Cauliflower Celery Chard Chickpea	Chicory Cucumber Fennel Leek Lentil Lettuce Marrow Onion Pea	Peppers (Mixed) Potato Quinoa Radish Rocket Shallot Soya Bean Spinach Squash	Sweet Potato Tomato Turnip Watercross Yuca
<b>Miscellaneous</b> 	Agar Agar Aloe Vera Cane Sugar	Carob Chestnut Cocoa Bean	Coffee Coia Nut Honey	Mushroom Tea (Black) Tea (Green)	Transglutaminase Yeast (Bakers) Yeast (Brewers)

Foods included in the FoodPrint® 200+ Panel

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